

Interton GN

PRO Series 
European Technology from Denmark 

Interton PRO Series Portfolio 2026

w.e.f 2nd February 2026



Works with
android 

Made for
iPhone | iPad | iPod

An uplifting sound experience

Users recommend Interton Presto for sound quality and hearing in noise. Start enjoying clear and easy conversations.

Styles designed to delight you

With the small, comfortable designs and our latest technology, you'll look forward to putting on your Interton Presto hearing aids.

Greater ways to stay connected

Connect with Bluetooth® Low Energy (LE) Audio and Auracast broadcast audio, the new global Bluetooth experience and assistive listening technology.



Hear your best anywhere

Focus on the details you enjoy or who you want to listen to while hearing clearly all around you.

Enjoy hands-free calls

Stream hands-free calls from iPhone, iPad and compatible Android™ mobile devices* directly to your Interton Presto hearing aids.

* Hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), iPad mini (6th generation), iPad (10th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later.
Bluetooth LE Audio works with Interton Presto and compatible mobile devices with Bluetooth 5.3 or later. See compatibility on interton.com/compatibility

What's your ideal style?

In the ear, behind the ear or custom made, there's a smaller-than-ever style that's great for you.

100% recyclable packaging

Your Interton Presto hearing aids come in packaging created from FSC certified paper and with 47% reduced CO2 emissions** as part of our focus on sustainability.

** Compared to previous hearing aid packaging from GN, CO2 reduction verified by Bureau Veritas.

Stylish wireless accessories

Designed for life's loud moments and made for Auracast™. Stream sound to boost your moments.

Entertainment and conversations

Hear TV clearly with the TV-Streamer+ and cut through the noise with the Multi-Mic+.

	Feature name	Presto Pro 9	Presto Pro 7	Presto Pro 5	Presto Pro 4
Audiological features	Number of channels	17	14	12	12
	M&RIE	✓	✓	✓	✓
	All access directionality ii	✓	✗	✗	✗
	Spatial sound	✓	✓	✓	✗
	Front focus	✓	✗	✗	✗
	Synchronised auto-steered directionality	✓	✓	✓	✓
	Environmental tuner ii	✓	✗	✗	✗
	Adaptive noise reduction	5 Settings (10 dB)	3 Settings (8 dB)	2 Settings (6 dB)	On / Off (3 dB)
	Microphone noise reduction	3 Settings	2 Settings	On / Off	On / Off
	Impulse noise reduction	3 Settings	3 Settings	On / Off	On / Off
	Adaptive wind noise reduction	3 Settings	2 Settings	On / Off	On / Off
	Frequency shifter	✓	✓	✓	✓
	Feedback manager plus (w/ music mode)	✓	✓	✓	✓
	Synchronised acclimatization manager	✓	✓	✓	✓
Functional features	Ear-to-ear communication	✓	✓	✓	✓
	Bluetooth® Auracast™	✓	✓	✓	✓
	Direct audio streaming (compatible ios** and android™ devices)	✓	✓	✓	✓
	Tv streamer 2, tv-streamer+, remote control, remote control 2, phone clip+	✓	✓	✓	✓
	Micro Mic, Multi Mic & Multi-Mic+	✓	✓	✓	✓
	Interton sound app	✓	✓	✓	✓
	Sound enhancer (interton sound app)	✓	✗	✗	✗
Online services	Interton remote fitting	✓	✓	✓	✓
	Interton remote fitting live	✓	✓	✓	✓
	Remote firmware updates	✓	✓	✓	✓
Fitting features	Fully flexible programs	4	4	4	4
	Pre-set feedback manager	✓	✓	✓	✓
	Datalogging	✓	✓	✓	✓
	Tinnitus sound generator	✓	✓	✓	✓
	Automatic receiver detection	✓	✓	✓	✓
	Wireless fitting with noahlink wireless	✓	✓	✓	✓



Sr. No.	Model name	Programming adapters	No. of channel	Battery types	Full on gain	Fitting range	MRP.in ₹
NEW Presto Pro Series 9 Made for iPhone iPad iPod Works with android Made for Android™ ASHA 							
Software - Interton fitting 2.0 & higher						Chipset - C 6	
1	Presto Pro 9 microRIE (2 Units + 1 Premium charger) MM / LP / MP / HP / UP 	Noahlink wireless	17	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	8,50,000 /-
2	Presto Pro 9 microRIE (1 Unit + 1 Standard charger) MM / LP / MP / HP / UP 	Noahlink wireless	17	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	4,35,000 /-
3	Presto Pro 9 HP BTE (2 Units + 1 Desktop charger) 	Noahlink wireless	17	Li-ion	80	30 - 110	8,50,000 /-
4	Presto Pro 9 HP BTE (1 Unit + 1 Desktop charger) 	Noahlink wireless	17	Li-ion	80	30 - 110	4,32,000 /-
5	Presto Pro 9 NP BTE (2 Units + 1 Desktop charger) 	Noahlink wireless	17	Li-ion	67	20 - 100	8,50,000 /-
6	Presto Pro 9 NP BTE (1 Unit + 1 Desktop charger) 	Noahlink wireless	17	Li-ion	67	20 - 100	4,32,000 /-
7	Presto Pro 9 ITC/ITE (2 Units + 1 Custom charger) LP / MP / HP / UP 	Noahlink wireless	17	Li-ion	50 / 60 / 70 / 78	20 - 110	8,50,000 /-
8	Presto Pro 9 ITC/ITE (1 Unit + 1 Custom charger) LP / MP / HP / UP 	Noahlink wireless	17	Li-ion	50 / 60 / 70 / 78	20 - 110	4,35,000 /-
9	Presto Pro 9 ITC LP / MP / HP / UP	Noahlink wireless	17	312	51 / 60 / 69 / 77	20 - 110	4,35,000 /-
10	Presto Pro 9 NWL CIC LP / MP / HP	Flex cable / HI PRO	17	10A	51 / 62 / 69	20 - 105	4,15,000 /-

Presto Pro Series 9 2 + 2 Years = 4 Years Total Warranty

NEW Presto Pro Series 7 Made for iPhone iPad iPod Works with android Made for Android™ ASHA 							
Software - Interton fitting 2.0 & higher						Chipset - C 6	
11	Presto Pro 7 microRIE (2 Units + 1 Premium charger) MM / LP / MP / HP / UP 	Noahlink wireless	14	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	4,60,000 /-
12	Presto Pro 7 microRIE (1 Unit + 1 Standard charger) MM / LP / MP / HP / UP 	Noahlink wireless	14	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	2,40,000 /-
13	Presto Pro 7 HP BTE (2 Units + 1 Desktop charger) 	Noahlink wireless	14	Li-ion	80	30 - 110	4,60,000 /-
14	Presto Pro 7 HP BTE (1 Unit + 1 Desktop charger) 	Noahlink wireless	14	Li-ion	80	30 - 110	2,38,000 /-
15	Presto Pro 7 NP BTE (2 Units + 1 Desktop charger) 	Noahlink wireless	14	Li-ion	67	20 - 100	4,60,000 /-
16	Presto Pro 7 NP BTE (1 Unit + 1 Desktop charger) 	Noahlink wireless	14	Li-ion	67	20 - 100	2,38,000 /-
17	Presto Pro 7 ITC/ITE (2 Units + 1 Custom charger) LP / MP / HP / UP 	Noahlink wireless	14	Li-ion	50 / 60 / 70 / 78	20 - 110	4,60,000 /-
18	Presto Pro 7 ITC/ITE (1 Unit + 1 Custom charger) LP / MP / HP / UP 	Noahlink wireless	14	Li-ion	50 / 60 / 70 / 78	20 - 110	2,40,000 /-
19	Presto Pro 7 ITC LP / MP / HP / UP	Noahlink wireless	14	312	51 / 60 / 69 / 77	20 - 110	2,20,000 /-
20	Presto Pro 7 NWL CIC LP / MP / HP	Flex cable / HI PRO	14	10A	51 / 62 / 69	20 - 105	2,20,000 /-

Presto Pro Series 7 2 + 2 Years = 4 Years Total Warranty

Sr. No.	Model name	Programming adapters	No. of channel	Battery types	Full on gain	Fitting range	MRP.in ₹
NEW Presto Pro Series 5 Made for iPhone iPad iPod Works with android Made for Android™ ASHA 							
Software - Interton fitting 2.0 & higher						Chipset - C 6	
21	Presto Pro 5 microRIE (2 Units + 1 Premium charger) MM / LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	2,70,000 /-
22	Presto Pro 5 microRIE (1 Unit + 1 Standard charger) MM / LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	1,50,000 /-
23	Presto Pro 5 HP BTE (2 Units + 1 Desktop charger) 	Noahlink wireless	12	Li-ion	80	30 - 110	2,60,000 /-
24	Presto Pro 5 HP BTE (1 Unit + 1 Desktop charger) 	Noahlink wireless	12	Li-ion	80	30 - 110	1,40,000 /-
25	Presto Pro 5 NP BTE (2 Units + 1 Desktop charger) 	Noahlink wireless	12	Li-ion	67	20 - 100	2,60,000 /-
26	Presto Pro 5 NP BTE (1 Unit + 1 Desktop charger) 	Noahlink wireless	12	Li-ion	67	20 - 100	1,40,000 /-
27	Presto Pro 5 ITC/ITE (2 Units + 1 Custom charger) LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	50 / 60 / 70 / 78	20 - 110	2,80,000 /-
28	Presto Pro 5 ITC/ITE (1 Unit + 1 Custom charger) LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	50 / 60 / 70 / 78	20 - 110	1,50,000 /-
29	Presto Pro 5 ITC LP / MP / HP / UP	Noahlink wireless	12	312	51 / 60 / 69 / 77	20 - 110	1,25,000 /-
30	Presto Pro 5 NWL CIC LP / MP / HP	Flex cable / HI PRO	12	10A	51 / 62 / 69	20 - 105	1,25,000 /-

Presto Pro Series 5 2 + 1 Years = 3 Years Total Warranty

NEW Presto Pro Series 4 Made for iPhone iPad iPod Works with android Made for Android™ ASHA 							
Software - Interton fitting 2.0 & higher						Chipset - C 6	
31	Presto Pro 4 microRIE (2 Units + 1 Premium charger) MM / LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	1,70,000 /-
32	Presto Pro 4 microRIE (1 Unit + 1 Standard charger) MM / LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	95,000 /-
33	Presto Pro 4 HP BTE (2 Units + 1 Desktop charger) 	Noahlink wireless	12	Li-ion	80	30 - 110	1,42,000 /-
34	Presto Pro 4 HP BTE (1 Unit + 1 Desktop charger) 	Noahlink wireless	12	Li-ion	80	30 - 110	80,000 /-
35	Presto Pro 4 NP BTE (2 Units + 1 Desktop charger) 	Noahlink wireless	12	Li-ion	67	20 - 100	1,42,000 /-
36	Presto Pro 4 NP BTE (1 Unit + 1 Desktop charger) 	Noahlink wireless	12	Li-ion	67	20 - 100	80,000 /-
37	Presto Pro 4 ITC/ITE (2 Units + 1 Custom charger) LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	50 / 60 / 70 / 78	20 - 110	1,70,000 /-
38	Presto Pro 4 ITC/ITE (1 Unit + 1 Custom charger) LP / MP / HP / UP 	Noahlink wireless	12	Li-ion	50 / 60 / 70 / 78	20 - 110	95,000 /-
39	Presto Pro 4 ITC LP / MP / HP / UP	Noahlink wireless	12	312	51 / 60 / 69 / 77	20 - 110	75,000 /-
40	Presto Pro 4 NWL CIC LP / MP / HP	Flex cable / HI PRO	12	10A	51 / 62 / 69	20 - 105	75,000 /-

Presto Pro Series 4 3 Years Total Warranty

Sr. No.	Model name	No. of channel	Battery types	Full on gain	Fitting range	MRP.in ₹		
NEW Presto Pro CROS Made for iPhone iPad iPod Works with android Made for Android™ ASHA 								
Software - Interton fitting 2.0 & higher		Programming Adapters - Noahlink wireless			Chipset - C 6			
41	Presto Pro 9 CROS microRIE Kit (1 Unit + 1 CROS* + 1 Premium charger)	MM / LP / MP / HP / UP		17	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	4,80,000 /-
42	Presto Pro 7 CROS microRIE Kit (1 Unit + 1 CROS* + 1 Premium charger)	MM / LP / MP / HP / UP		14	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	2,85,000 /-
43	Presto Pro 5 CROS microRIE Kit (1 Unit + 1 CROS* + 1 Premium charger)	MM / LP / MP / HP / UP		12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	1,90,000 /-
44	Presto Pro 4 CROS microRIE Kit (1 Unit + 1 CROS* + 1 Premium charger)	MM / LP / MP / HP / UP		12	Li-ion	62 / 65 / 70 / 73 / 80	20 - 110	1,40,000 /-
45	Presto Pro CROS microRIE (1 CROS" Transmitter)	LP / MP / HP			Li-ion			59,300 /-
46	Presto Pro CROS NP BTE Transmitter (End use will pair it with a Presto 88 or 77 BTE monaural kit)				Li-ion			59,300 /-

Presto Pro CROS

1 Year Warranty on Transmitter microRIE



	Feature Name	Move Pro 6	Move Pro 4	Move Pro 3	Move Pro 2	Move Pro 1
Sound quality	Warp compression (WDRC), number of channels	17	12	8	6	4
	Adaptive noise reduction	✓	✓	✓	✓	✗
	Adaptive wind noise reduction	✓	✓	✓	✓	✗
	Impulse noise reduction	✓	✓	✗	✗	✗
	Microphone noise reduction	✓	✓	✓	✓	✓
	Environmental gain tuner	✓	✗	✗	✗	✗
	Environmental classifier	✓	✓	✓	✓	✓
Speech understanding	Integrated directionality	✓	✗	✗	✗	✗
	Automatic beamwidth	✓	✗	✗	✗	✗
	Combined directionality	✓	✓	✗	✗	✗
	Synchronised auto-steered directionality	✓	✓	✗	✗	✗
	Selectable beamwidth	✓	✓	✓	✗	✗
	Auto-steered directionality	✓	✓	✓	✓	✗
	Speech-focused directionality	✓	✓	✓	✓	Omni / Fixed
Feedback management	Feedback manager plus	✓	✓	✓	✓	✗
	Music mode	✓	✓	✗	✗	✗
	Preset feedback manager	✓	✓	✓	✓	✓
Apple & android connectivity / rechargeability	Direct audio streaming (MFI, android™)	✓	✓	✓	✗	✗
	Li-ion rechargeable battery	✓	✓	✓	✗	✗
Onboarding	Synchronised acclimatization manager	✓	✓	✓	✗	✗
	Acclimatization manager	✓	✓	✓	✓	✗
Convenience features	Ear to ear communication (push button, volume control)	✓	✓	✓	✗	✗
	Power-on-delay	✓	✓	✓	✓	✗
	Autophone	✓	✓	✓	✓	✓
Wireless accessories & apps	TV-Streamer 2	✓	✓	✓	✓	✗
	Remote control 2	✓	✓	✓	✓	✗
	Phone clip 2	✓	✓	✓	✓	✗
	Multi Mic	✓	✓	✓	✓	✗
	Interton Sound™ app	✓	✓	✓	✓	✓
	Remote firmware update	✓	✓	✓	✓	✓
Tinnitus solutions & Catalog	Tinnitus sound generator	✓	✓	✓	✓	✗
	Datalogging	✓	✓	✓	✓	✓

Receiver-In-the-Ear (RIE) hearing aids



MicroRIE



Rechargeable RIE



Invisible-In-the-Canal (IIC)



Completely-In-the-Canal (CIC)



Completely-In-the-Canal (CIC-W)



In-the-Canal (ITC)

Custom-made hearing aids

Behind-the-Ear (BTE) Hearing aids



SP BTE



Power BTE

Sr. No.	Model name	Programming adapters	No. of channel	Battery types	Full on gain	Fitting range	MRP.in ₹
Interton Move Pro 6				Made for iPhone iPad iPod	Works with android	Made for Android™ ASHA	
Software - Interton fitting 2.0 & higher				Chipset - C 5 (Coyote)			
1	Move Pro 661 DRWC RIE (2 Units + 1 Premium Charger) LP / MP / HP / UP 	Noahlink wireless	17	Li-Ion	62 / 67 / 74 / 82	20 - 110	1,40,000 /-
2	Move Pro 661 DRWC RIE (1 Unit + 1 Desktop Charger) LP / MP / HP / UP 	Noahlink wireless	17	Li-Ion	62 / 67 / 74 / 82	20 - 110	74,000 /-
3	Move Pro 690 DW SP BTE	Noahlink wireless	17	675	86	60 - 120	80,000 /-
4	Move Pro 680 DWH HP BTE Metal Hook	Noahlink wireless	17	13	80	30 - 110	60,000 /-
5	Move Pro 670 DW BTE Thin tube / Closed NP BTE	Noahlink wireless	17	13	63 / 66	20 - 100	60,000 /-
6	Move Pro 6 DW ITC MP / HP / UP	Noahlink wireless	17	312	59 / 69 / 78	20 - 110	60,000 /-
7	Move Pro 6 NWL CIC LP / MP / HP	HI PRO	17	10A	49 / 59 / 69	20 - 105	60,000 /-
Interton Move Pro 4				Made for iPhone iPad iPod	Works with android	Made for Android™ ASHA	
Software - Interton fitting 2.0 & higher				Chipset - C 5 (Coyote)			
8	Move Pro 461 DRWC RIE (2 Units + 1 Premium charger) LP / MP / HP / UP 	Noahlink wireless	12	Li-Ion	62 / 67 / 74 / 82	20 - 110	1,12,000 /-
9	Move Pro 461 DRWC RIE (1 Unit + 1 Desktop charger) LP / MP / HP / UP 	Noahlink wireless	12	Li-Ion	62 / 67 / 74 / 82	20 - 110	61,000 /-
10	Move Pro 490 DW SP BTE	Noahlink wireless	12	675	86	60 - 120	60,000 /-
11	Move Pro 461 DRW RIE LP / MP / HP / UP	Noahlink wireless	12	312	62 / 67 / 74 / 82	20 - 110	42,000 /-
12	Move Pro 462 DRW RIE LP / MP / HP / UP	Noahlink wireless	12	13	62 / 67 / 74 / 82	20 - 110	42,000 /-
13	Move Pro 480 DWH HP BTE Metal Hook	Noahlink wireless	12	13	80	30 - 110	42,000 /-
14	Move Pro 470 DW BTE Thin tube / Closed NP BTE	Noahlink wireless	12	13	63 / 66	20 - 100	42,000 /-
15	Move Pro 4 DW ITC MP / HP / UP	Noahlink wireless	12	312	59 / 69 / 78	20 - 110	50,000 /-
16	Move Pro 4 NWL CIC LP / MP / HP	HI PRO	12	10A	49 / 59 / 69	20 - 105	50,000 /-



Rechargeable
Hearing Aids



Apple & Android
ASHA Connectivity



Water resistant
Isolate nano-coating



Interton
Sound app


Sr. No.	Model name		Programming adapters	No. of channel	Battery types	Full on gain	Fitting range	MRP.in ₹	
Interton Move Pro 3					Made for iPhone iPad iPod		Works with android	Made for Android™ ASHA	
Software - Interton Fitting 2.0 & Higher					Chipset - C 5 (Coyote)				
17	Move Pro 361 DRWC RIE (2 Units + 1 Desktop charger)	LP / MP / HP / UP 	Noahlink wireless	8	Li-Ion	62 / 67 / 74 / 82	20 - 110	92,000 /-	
18	Move Pro 361 DRWC RIE (1 Unit + 1 Desktop charger)	LP / MP / HP / UP 	Noahlink wireless	8	Li-Ion	62 / 67 / 74 / 82	20 - 110	52,000 /-	
19	Move Pro 390 DW SP BTE		Noahlink wireless	8	675	86	60 - 120	41,000 /-	
20	Move Pro 361 DRW RIE	LP / MP / HP / UP	Noahlink wireless	8	312	62 / 67 / 74 / 82	20 - 110	31,000 /-	
21	Move Pro 362 DRW RIE	LP / MP / HP / UP	Noahlink wireless	8	13	62 / 67 / 74 / 82	20 - 110	31,000 /-	
22	Move Pro 380 DWH HP BTE Metal hook		Noahlink wireless	8	13	80	30 - 110	31,000 /-	
23	Move Pro 370 DW BTE Thin tube / closed NP BTE		Noahlink wireless	8	13	63 / 66	20 - 100	31,000 /-	
24	Move Pro 3 DW ITC	MP / HP / UP	Noahlink wireless	8	312	59 / 69 / 78	20 - 110	32,000 /-	
25	Move Pro 3 NWL CIC	LP / MP / HP	HI PRO	8	10A	49 / 59 / 69	20 - 105	32,000 /-	
Interton Move Pro 2									
Software - Interton Fitting 2.0 & Higher					Chipset - C 5 (Coyote)				
26	Move Pro 290 DW SP BTE		Noahlink wireless	6	675	86	60 - 120	30,700 /-	
27	Move Pro 261 DRW RIE	LP / MP / HP / UP	Noahlink wireless	6	312	62 / 67 / 74 / 82	20 - 110	26,000 /-	
28	Move Pro 262 DRW RIE	LP / MP / HP / UP	Noahlink wireless	6	13	62 / 67 / 74 / 82	20 - 110	26,000 /-	
29	Move Pro 280 DWH HP BTE Metal Hook		Noahlink wireless	6	13	80	30 - 110	28,500 /-	
Interton Move Pro 1									
Software - Interton Fitting 2.0 & Higher					Chipset - C 5 (Coyote)				
30	Move Pro 180 DWH HP BTE Metal Hook		Noahlink wireless	4	13	80	30 - 110	23,000 /-	

Note : D - Directional | R - RIE | W - Wireless | C - Chargeable

3 Years Warranty Applicable to Hearing Aid Models for which the warranty is not specified.



Interton Features Overview




Integrated Directionality

Adaptive microphone feature used in binaural fittings. Continuous ear-to-ear exchange of data about the environment optimises microphone settings in each ear for a wide variety of everyday listening situations.

Benefit

Users can focus on the sounds that are important to them without losing awareness of what's going on around them, so they can effortlessly follow conversations and shift attention naturally in any sound environment, all without the need for any manual adjustments.




Adaptive Noise Reduction

A unique and fast-acting spectral subtraction technology that precisely identifies speech and noise and reduces unwanted noise without degrading the speech signal.

Benefit

Noise reduction will provide comfort in noisy listening situations without compromising speech understanding.




Auto-Steered Directionality

Hearing aids analyse the sound environment and switch microphone settings between omnidirectional and directional. In directional mode, the beam width can adapt to maximise listening in noise. The switching happens so smoothly that users don't notice the change.

Benefit

Users can benefit from automatic switching for better speech understanding when needed without any need to switch programs manually. The transitions between microphone modes are performed smoothly to provide a more natural sound experience.




Speech-Focused Directionality

An adaptive directional feature, which changes directional polar patterns based on the surrounding noise source for the optimal directional benefit.

Benefit

Users can understand speech in front of them better even in dynamic sound environments.




Feedback Manager Plus

A feedback cancellation feature that uses two cancellation filters to eliminate feedback without reducing amplification or compromising sound quality. Feedback Manager Plus also has the ability to eliminate feedback from sudden movements or objects placed near the ear.

Benefit

Best-in-class feedback cancellation enables users to enjoy the sound quality without experiencing feedback or distorted sound in any situation.



Environmental Gain Tuner

Environmental Gain Tuner optimises gain settings for better listening experiences in five different environments. The hearing aids analyse the surrounding environment and automatically apply preset gain settings as the user moves from one acoustic environment to another. The gain settings can be customised for each environment based on personal preference.

Benefit

Users can benefit from automatic adjustments in gain without the need to make any adjustments in volume themselves. It's the perfect feature for active people who encounter a wide variety of listening situations. The customisable gain settings inside the feature ensure that each user experiences the best settings specific to them and the situations they encounter.




Acclimatization Manager

For some users, adapting to prescribed gain can be a challenge, especially when the users are new to hearing aids. Acclimatization Manager gradually increases gain over time to eventually achieve prescribed gain. The starting point can be set as low as 60% of the desired end gains, and the amount of time over which the gain increases can be customised to each user.

Benefit

Acclimatization Manager allows new users to adapt to amplification comfortably and seamlessly. It can also be beneficial in cases where a new user is unable to return to the dispenser for follow-up for an extended period of time.



Combined Directionality

Microphone feature used in binaural fittings. Microphone settings are set at the time of hearing aid programming by hearing care professional in each ear for a wide variety of everyday listening situations.

Benefit

Users will be able to focus on speech without losing the awareness of environment, more like a nature intended hearing.

Cable & fitting guide


	 (NOAHLINK WIRELESS)	 (HI-PRO CABLE/CS44) (ADAPTER, FLEX, 3 PIN, CS53)	 (HI-PRO CABLE/CS44) (PROGRAMMING ADAPTOR 2)	 IIC Programming	 Fitting Software
Interton Presto Pro					Interton Fitting 2.0 or Higher
Micro RIE	●				
BTE	●				
ITC / ITE	●				
CIC		●			
IIC				●	
Interton Move Pro					Interton Fitting 2.0 or Higher
RIE Rechargeable	●				
RIE 312 Battery	●				
RIE 13 Battery	●				
Normal Power BTE (70)	●				
High Power BTE (80)	●				
Super Power BTE (90)	●				
Custom WL	●				

How to pair your hearing aids with your Android™ device



1. Turn your hearing aids off.
2. Go to google play on your phone search for interton sound app.
3. Download interton sound app & install.
4. Go the settings menu and check that bluetooth is turned on.
follow instructions on the screen.
5. Turn your hearing aids on.
6. Click to open the interton sound app follow instructions
on the screens to connect.

How to pair your hearing aids with your iOS device



1. Turn your hearing aids off.
2. Turn on bluetooth by going to the settings menu on your iOS device.
3. Go back to the main settings menu and then, tap on general.
4. In the general menu, tap on hearing aid.
5. Turn your hearing aids on.
6. Turn your hearing aids on your iphone / ipad / ipod will detect your hearing aids.
7. Tap on pair.
8. Go to the app store, install interton sound app, follow instructions to connect to the app.

Charger feature	Premium charger (C-1)	Standard charger (C-2)	Custom charger (C-4)	Desktop charger (C-3)
Onboard power bank	YES	NO	NO	NO
LED	YES	YES limited functionality	YES limited functionality	YES limited functionality
Lid	YES	YES	YES	NO
Removable cable	USB-C	USB-C	USB-C	USB-C
Inserts (CFF)	CFF 08 for Presto MicroRIE	CFF 08 for Presto MicroRIE	Custom	CFF70 or CFF80

Premium charger



Premium Charger:

- Wireless charging solution
- Acts as powerbank
- Once charged will last for 10 Days
- Charges hearing aids for 3 times
- Fast charging enabled
- Easy to use LED indicators for battery life indication

Standard charger



Standard Charger:

- wireless charging solution
- once charged will last for 30 hours
- fast charging enabled
- comes with a lid to carry the charger

Custom charger



Custom Charger:

- Wireless charging solution
- Once charged will last for 24 hours
- Fast charging enabled
- Comes with a lid to carry the device
- The industry's only customised charger to ensure reliable charging

Desktop charger



Desktop Charger:

- Wireless charging solution
- Once charged will last 30 hrs
- Fast charging enabled

Chargers

MRP. in ₹

1	Premium Charger (Presto Pro MicroRIE, Move Pro - RIE 4 & 6)	15,500 /-
2	Standard Charger (Presto Pro MicroRIE)	10,500 /-
3	Custom Charger (Presto Pro Rechargeable ITC)	10,500 /-
4	Desktop Charger (Presto Pro BTE, Move Pro - RIE 3, 4 & 6)	8,500 /-

1 Year Warranty on all types of Chargers



TV-Streamer+

Multi Mic

TV-Streamer 2

Multi-Mic+

Wireless accessories

MRP. in ₹

1	TV-Streamer+ / TV-Streamer 2	17,700 /-
2	Multi-Mic+ / Multi Mic	17,700 /-

1 Year Warranty on Wireless Accessories

Terms & Conditions

1. Interton PRO hearing aids carry 2+1 Years Warranty on MOVE Pro 1,2,3,4,6 series product & Presto Pro 4 & 5 series products.
2. Interton Presto Pro 7 & Presto Pro 9 series Hearing aids will carry 2+2 Years Warranty.
3. Warranty will be considered only against any manufacturing defects. Warranty will not be honoured if hearing aids are damaged due to mishandling, physical damage, battery leakage, earwax, liquids, oils, etc.
4. Accessories (Chargers, RIE Receiver's, Remote control, TV-streamers, Multi Mic) carry a 1-year warranty against manufacturing defects. Warranty will not be honoured if accessories are damaged due to mishandling, physical damage, battery leakage, earwax, liquids, oil, etc.
5. Power adaptor for some accessories will have to be procured separately.
6. Custom hearing aids prices include shelling charges.
7. HP/UP RIE Receiver's customization charges & RIE mould charges are not included in the RIE pricing. It will be available at an additional cost.
8. RIE / BTE instruments will be supplied in standard beige colour. Other colour housing options for RIE / BTE models are subject to a delivery lead time of 3 weeks.
9. Prices of the colour conversion kit may vary depending on hearing aid model.
10. Hearing aids are exempted from GST.
11. GST is applicable on accessories. Retail prices of accessories are inclusive of GST.
12. The pricing and product specifications are subject to change without prior notice.
13. This price list supersedes all previous price lists.
14. We recommend all our valued customers to seek a hearing care professional's opinion on the best-suited hearing instruments. Online purchase may not necessarily give you the desired results.
15. Warranty card must accompany all warranty repair cases, repair cases received without warranty card will be treated as out of warranty.
16. This Price list is effective from 2nd February 2026.

Find us on :



facebook.com/Arphielectronics



twitter/Arphielectronics



instagram.com/Arphielectronics



linkedin/Arphielectronics



<https://www.youtube.com/channel/UCSIN0p3-4KkaQXak7rCZZpw>

Discover your perfect hearing solution today



Interton GN
.....

Website: www.interton.com

Marketed by:

Arphi® ELECTRONICS PVT. LTD.

BRANCHES : Mumbai • New Delhi • Chennai • Kolkata

Website: www.arphi.com